


Increased Functionality of Optical Fibers for Life-Science Applications

Aziza Sudirman
 Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH



Motivation and Aim

Question:

- Why use optical fibers?
- How can life-science applications benefit from optical fibers?

Answer:

- Ⓢ Excellent waveguides.
- Ⓢ Micrometer-size cross section.
- Ⓢ Thin, flexible and can be made very long.
- Ⓢ However, cannot offer same complexity as lab-on-a-chip.

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Motivation and Aim

But, optical fibers offer advantages such as...

- Ⓢ Allow for real-time measurements.
- Ⓢ Can be inserted into the body.
- Ⓢ Core-size allows for light delivery to single biological species.

The aim is to work towards optical analysis of biological species at a single-cell level for in-vivo studies.

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH


Outline

- Four different areas (proof-of-principle studies):
 1. **Combining** laser light delivery with reflectometry and fluid collection in an optical fiber.
 2. **Developing** all-fiber optofluidic components.
 3. **Reducing** background luminescence in fiber-based fluorescence sensing.
 4. **Detecting & collecting** micrometer-size particles using a microstructured fiber.
- Summary and final remarks.

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Introduction



<p>Microstructured fiber:</p> <ul style="list-style-type: none"> • 125 μm outer-diameter. • 8 μm core. • Longitudinal holes. 	<p>Capillary:</p> <ul style="list-style-type: none"> • 125 μm outer-diameter. • Standard telecom fiber. • Various inner-diameter. 	<p>Carbon-coated fiber:</p> <ul style="list-style-type: none"> • 125 μm outer-diameter. • Standard telecom fiber. • ~20 nm thick layer of carbon.
--	---	---

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Four different areas



Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

What do we want to solve?

- 1) Method for positioning of fiber-tip.
Solution: Reflectometry.
- 2) Removal of sample that is "in the way".
Solution: Delivery of high-power laser light for ablation using optical fiber.
- 3) Perform biopsy - collecting fluid sample.
Solution: Exploit the holes of microstructured fiber.

Laser ablation, reflectometry and microfluidics

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Laser ablation, reflectometry and microfluidics

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Four different areas

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

What do we want to solve?

When dealing with coupling light and fluid in a fiber...

- 1) How to handle fluid in a controlled way without disturbing light guidance.
- 2) Eliminate issues with liquid evaporation/dripping.
- 3) Prevent forming of liquid meniscus:

Solution: All-in-fiber components that are user-friendly.

Optofluidic components: "Billys"

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Optofluidic components: "Billys"

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Optofluidic components

Allows for interaction of light and fluid...

Small inner-diameter capillary allows for:

- ☺ Stronger coupling between light and fluid.
- ☺ Low-loss light coupling to core of microstructured fiber.

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Optofluidic components

All-fiber component:

- ☺ Fluid-filled or continuous flow of fluid.
- ☺ Reflected or transmitted signal for measurements.
- ☺ Entirely sealed, with no risks for fluid contamination.

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Optofluidic components

Replace... with...

Traditional microfluidic component Fiber-based component

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Four different areas

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

What do we want to solve?

When dealing with low-level signal sensing.

1) Luminescence from fiber is a source for noise.
Solution: Fiber with carbon-coating.

Reducing fiber luminescence

Characterizing luminescence from four different fiber types

Compared four fiber types:

- Standard telecom fiber with Ge-doped core
- Pure silica fiber
- Carbon-coated standard telecom fiber with Ge-doped core
- Carbon-coated pure silica fiber

depr. cladi.

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Reducing fiber luminescence

Detecting Raman signal from organic solvents

Standard telecom fiber Carbon-coated fiber

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Four different areas

Detecting & collecting
Fluorescent particles

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

What do we want to solve?

- 1) Detect single-particles of interest.
Solution: Laser-induced fluorescence sensing.
- 2) Isolate detected particles.
Solution: Exploit holes of microstructured fiber.

Detecting and collecting fluorescent particles

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Detecting and collecting fluorescent particles

Fluorescence intensity (arb. units)

Time (sec)

Signal from fluorescent bead 15 μm away from fiber-tip

Collected bead inside fiber hole

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Summary

- Laser light delivery for ablation combined with positioning of fiber-tip based on reflectometry.
- Collection of fluid and particles into the hole(s) of microstructured fiber.
- All-fiber components for combining light and fluid.
- Carbon-coated fiber for reducing luminescence from polymer coating of fiber.
- Reflected signal for detecting of Raman and fluorescence.

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH

Final remarks



Thank you for your attention!

Sudirman A., Doctoral Thesis Defence, May 21, 2014 in FDS, Albanova KTH